

# ABON™ Tumor Marker Rapid Tests

## Early Diagnosis of Dreadful Tumors

- STAT Diagnosis
- Simple Procedure
- Accurate Results



# ABON™ Tumor Marker Rapid Tests

## Early Diagnosis of Dreadful Tumors

### STAT Diagnosis

- Results in minutes
- Facilitates treatment decisions early

### Simple Procedure

- Easy to perform and interpret
- No additional equipment require


### Accurate Results

- High sensitivity and specificity
- Built-in procedural control

## Test Procedure


### CEA Whole Blood Strip Test

- Add specimen and buffer
- Read results in minutes




### AFP CEA Whole Blood Device Test

- Add specimen and buffer
- Read results in minutes



### FOB Strip/Device Test

- Collect specimen with collection stick
- Agitate vigorously to mix specimen and buffer
- Immerse or add extracted specimen to the device and read results at 5 minutes



### PSA Whole Blood Strip/Device Test

- Add specimen and buffer
- Read results at 5 minutes

Semi-Quantitative




Qualitative



## Ordering information

Product Name	Specimen	Format	Catalog No.	Read Time	Sensitivity	Tests per kit
AFP Alpha-Fetoprotein	Whole Blood/Serum/Plasma	Device	TAF-402 ✓	10 min.	10 ng/mL	40
CEA Carcinoembryonic Antigen	Whole Blood/Serum/Plasma	Strip	TCE-401	5 min.	5ng/mL	50
		Device	TCE-402 ✓			40
FOB Fecal Occult Blood	Feces	Strip	TFO-601 ✓†	5 min.	50 ng/mL 6 µg/g	25
		Device	TFO-601S ✓			40
		Device	TFO-602 ✓			25
		Strip	TPS-401			50
PSA Prostate Specific Antigen Semi-Quantitative	Whole Blood/Serum/Plasma	Device	TPS-402 ✓	5 min.	4 ng/mL (R=10 ng/mL) 3 ng/mL (R=10 ng/mL)	40/25
			TPS-U402 ✓			
PSA Prostate Specific Antigen Qualitative	Whole Blood/Serum/Plasma	Device	TPS-UB402	5 min.	3 ng/mL	40

† US FDA Cleared ✓ CE Marked CE



**Inverness Medical Innovations Hong Kong Limited**  
 #15th Floor, Building B, Modern International Plaza,  
 159 Tianmu Shan Road, Hangzhou 310007, P.R.China  
 Please visit our website for details: [www.imihongkong.com](http://www.imihongkong.com)  
 E-mail: [IMHK.info@invmed.com](mailto:IMHK.info@invmed.com)

**TECNIGEN S.A.**  
 Av. Zañartu 1060  
 Ñuñoa, Santiago  
 Chile  
[www.tecnigen.cl](http://www.tecnigen.cl)  
[info@tecnigen.cl](mailto:info@tecnigen.cl)

